


## Disposable Ligation Device

# Disposable Ligation Device

---

## Intended Use

The product is used in combination with the endoscope to release the ligation loops for controlling or preventing sudden bleeding after excising pedunculated polyps.



Multiple choices of loop number

Rotatable loop hook



## Features

- Various nylon loop sizes designed for different focus
- Separated structure of releasing device and loops
- Flexible loop number combination
- Rotatable loop hook
- Self-lubricating sheath



Beijing ZKSK Technology Co., Ltd.

Building 9, 6 & No.6 Yuan Hengye North 7th street, Yongle Economic Development Zone, Tongzhou District, Beijing 101105, China

Tel: 86 10 8052 5660 Email: shawna@bjzksk.com

Website: [www.bjzksk.com/en](http://www.bjzksk.com/en) <http://bjzksk.en.alibaba.com>

Model	Opening diameter of ligation ring(mm) ±5 m m	Conveyor working length(c m) ±50 m m	Number of ligation rings (pcs/set)	Lead-in sheath O.D(mm) ± 0.15mm	Casing O.D (m m) ±0.3m m	Effective length of ligation ring(m m) ±5 m m
HE15-165-1	15	165	1	2.55	1.9	35
HE15-165-2		165	2	2.55	1.9	35
HE15-165-3		165	3	2.55	1.9	35
HE15-195-1		195	1	2.55	1.9	35
HE15-195-2		195	2	2.55	1.9	35
HE15-195-3		195	3	2.55	1.9	35
HE15-230-1		230	1	2.55	1.9	35
HE15-230-2		230	2	2.55	1.9	35
HE15-230-3		230	3	2.55	1.9	35
HE20-165-1	20	165	1	2.55	1.9	40
HE20-165-2		165	2	2.55	1.9	40
HE20-165-3		165	3	2.55	1.9	40
HE20-195-1		195	1	2.55	1.9	40
HE20-195-2		195	2	2.55	1.9	40
HE20-195-3		195	3	2.55	1.9	40
HE20-230-1		230	1	2.55	1.9	40
HE20-230-2		230	2	2.55	1.9	40
HE20-230-3		230	3	2.55	1.9	40
HE20-165-1	30	165	1	2.55	1.9	40
HE30-165-1		165	1	2.55	1.9	55
HE30-L165-2		165	2	2.55	1.9	55
HE30-L165-3		165	3	2.55	1.9	55
HE30-L195-1		195	1	2.55	1.9	55
HE30-L195-2		195	2	2.55	1.9	55
HE30-L195-3		195	3	2.55	1.9	55
HE30-L230-1		230	1	2.55	1.9	55
HE30-L230-2		230	2	2.55	1.9	55
HE30-L230-3	230	3	2.55	1.9	55	
HE40-L165-1	40	165	1	2.55	1.9	72
HE40-L165-2		165	2	2.55	1.9	72
HE40-L165-3		165	3	2.55	1.9	72
HE40-L195-1		195	1	2.55	1.9	72
HE40-L195-2		195	2	2.55	1.9	72
HE40-L195-3		195	3	2.55	1.9	72
HE40-L230-1		230	1	2.55	1.9	72
HE40-L230-2		230	2	2.55	1.9	72
HE40-L230-3		230	3	2.55	1.9	72